

# Release Notes

## StemLab 3.5.4

**STEMSOFT is very pleased to announce the release of StemLab 3.5.4**

### **IMPORTANT NOTES**

For those customers using BMTserve, please note that this StemLab upgrade is only compatible with **BMTserve 3.2.1 2008080808**

If you have any difficulties or questions, please do not hesitate to contact us. Call the STEMSOFT Support Line at 1-800-671-3234 or email [support@STEMSOFT.com](mailto:support@STEMSOFT.com).

### **TESTED CONFIGURATIONS**

- Microsoft SQL Server 2000 SP4 and 2005 SP3
- Oracle 9i and 10g
- Microsoft Windows XP Service Pack 3

### **ENHANCEMENTS**

StemLab now supports SOP paths to be URLs

- All existing rules of how StemLab handles paths for files are still in place.
- A URL must begin with either “http” or “https” and can point to a web page (e.g. <http://www.google.com>) or to a file (e.g. <http://yourservername/SOP.doc>)
- The URL will open in the default web browser

StemLab 3.5.4 also features several improvements to the Reporting Views. Specific issues addressed are documented below. Please see the Reporting View Guide for documentation of the reporting views. If you are currently using a version older than 3.5.3 please refer to the previous release notes for additional new features, including ISBT128 compliant functionality.

StemLab™ 3.5.4

# Release Notes

## StemLab 3.5.4

### ISSUES RESOLVED IN THIS RELEASE

#### **Issue 1967**

Previous Count for Pooled Product not calculated incorrectly when enumerations Assumed Value is set to "COUNT"  
Enumeration calculation logic re-written to address issue and improve the performance of enumeration calculations.

#### **Issue 2092**

Use a division process (auto-div or manual div) and divide a product into 10 or more fractions. Add inventory info for each fraction, and select Pool. The Pool window will display the number of fractions minus 1, with the first fraction missing.

#### **Issue 2098**

Slow performance when loading the Processing Step details.  
When loading the processing step details either in the process record or in maintenance mode, the grid will load line by line and can take.

#### **Issue 2007**

Runtime Error when changing Product Release State in a Work list  
There were ported cases of a runtime error when user changed the Product Release State of a work list.

#### **Issue 2090**

Performance improvements made to accessioning and process record tree

#### **Issue 1987**

Fraction Title Inheritance Not Functioning when segments are used.  
Title inheritance works as specified when stand-alone process is applied; does not work properly when segment is applied

#### **Issue 2140**

Resolved timeout issues during product summary load

### REPORTING VIEWS:

# Release Notes

## StemLab 3.5.4

**Issue 2171**

V\_R\_PatLabTest view does not populate with data.  
Issue affected only Oracle.

**Issue 2109**

In some cases the V\_R\_SL\_FREEZER\_LOCATION view did not match with actual location of cryopreserved products

**Issue 2100**

New column added to V\_R\_SL\_FREEZER\_LOCATION reporting view  
*FRACTION\_TITLE*

**Issue 2112**

New column added to V\_R\_SL\_SOURCE\_PRODUCTS  
*IDPROJECT*

Original column IDPROJECT\_COLLECTION left for backward compatibility

**Issue 2001**

V\_R\_SL\_ISBT128 needs Collection Facility

The ISBT label standard requires that you print the collection facility name and address. This was not included in the current view. Additional columns were added

*COLLECTION\_FACILITY\_CITY*  
*COLLECTION\_FACILITY\_COUNTRY*  
*COLLECTION\_FACILITY\_NAME*  
*COLLECTION\_FACILITY\_PROVINCE*

**Issue 1993**

New Columns added to the V\_R\_SL\_SOURCE\_PRODUCTS reporting view

*DONOR\_WEIGHT; MOBILIZATION\_TYPE*  
*MOBILIZATION\_DATE*

**Issue 1992**

New column added to the V\_R\_SL\_FINALPRODUCT\_TERMINAL reporting view

*PROCESSING\_DATE*